

5

reasons

you need

IceFORCE™ 2.1 CX L and IcePearl™ 2.1 CX L

Cryoablation Needles

in your ablation toolkit

Superior Freezing Performance

Patented, enhanced Joule-Thomson technology drives a 2.1mm shaft without compromising 2.4 mm iceball size or volume



Economic, Eco-efficient Operation

A categorical improvement over helium-thaw needles, 2.1 CX technology pairs advanced freezing capability with helium-free thawing technology to reduce total gas cost by up to 65%



Exclusive Thaw Alternatives

Proprietary helium-free thaw options defy conventional thawing techniques allowing up to 18% faster thaw time and quick needle withdrawal



Unique Cautery Functionality

Innovative cautery for track ablation provides confidence by creating a small zone of cellular necrosis around the needle track



Optimized Needle Design

A strong yet light-weight needle design takes weight out of the needle assembly to prevent torque and provide unmatched needle placement flexibility

